


## UCderm impression 1 x 1 ml

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**EN**

**Description:** UCderm impression is a sterile, isotonic, pyrogen-free, colourless and transparent gel made of highly purified, crosslinked hyaluronic acid of non-animal origin in a physiological phospat buffer. UCderm impression contains 0.3% lidocaine hydrochloride.

**Content:** Each package contains one single-use glass syringe, pre-filled with 1 ml UCderm impression and sterilised by moist heat. In addition, there is one single-use 30G sterile needle to be used only for injecting UCderm impression, an instruction leaflet, an injection instruction, a patient card and two labels in order to ensure traceability in it.

**Compositions:** UCderm impression contains:

- 26 mg/ml crosslinked Hyaluronic acid
- 3 mg/ml Lidocaine hydrochloride
- Phosphate buffer

The amount of gel in each syringe is 1 ml.

**Indications/Intended use:** UCderm impression is a sterile, injectable, absorbable implant for hydrating soft tissues in the face, neck and décolleté. This implant is indicated to balance the skin’s natural age-related loss of hyaluronic acid, to refine skin structure and to reduce fine lines caused by dryness. UCderm impression is intended for the planar treatment of skin areas. The product is intended to be used via upper dermis injection described in the enclosed injection instruction. The lidocaine contained in the product is meant to reduce the patient’s local pain during and for a short time after the injection.

**Contraindications:** Do not apply UCderm impression in the following cases:

- In youths under 18 years of age.
- During pregnancy and breast-feeding.
- In case of incipient or existing influenza infections.
- In case of acute systemic infections.
- In case of skin inflammations (e.g. acne, herpes, psoriasis)
- In case of a known hypersensitivity to sodium hyaluronate or lidocaine hydrochloride or amide-based local anaesthetics.
- In case of known autoimmune diseases.
- In patients with blood clotting disorders or patients on anti-coagulation medication (e.g. aspirin, non-steroidal anti-inflammatory drugs).
- In case of acute articular rheumatism with cardiac localization.
- In patients with a tendency to hypertrophic scarring
- In case of untreated epilepsy.
- In patients who suffer from porphyria.
- During treatment with drugs tending to slow down or inhibit the liver metabolism (e.g. cimetidine, beta-blockers).
- In case of cardiac arrhythmia.

**Precautions for use:** UCderm impression must not be injected in the glabellar region (forehead) and periorbital area (eyelids, eye shadows, crow’s feet). The application of UCderm impression for reducing the tear trough in the under-eye area is an exception. The application has to be performed only by specialists particularly trained in this technique who have a precise knowledge of the physiology of this certain area. For this area, the use of a blunt cannula is recommended.

When treating facial areas with limited soft tissue support or soft tissue coverage, such as the periorbital area, special care should be taken to avoid the formation of palpable lumps, discoloration or granulomas.

UCderm impression must not be injected into tendons, muscles, bones, joints or eyes.

UCderm impression must not be injected into blood vessels. Introduction of product into the vasculature may lead to embolization or arterial ischemia or infarction. Therefore, the product should always be injected slowly and at the lowest possible pressure.

UCderm impression must not be used simultaneously with other aesthetic medicine techniques such as laser treatments, chemical peelings or dermabrasion. If inflamed skin areas occur after the treatment with these techniques, injection into those areas is not recommended. The treated skin area has to heal completely first.

UCderm impression must not be injected into skin areas which had been subjected to a treatment with permanent implants, implants of animal origin or implants containing other substances than hyaluronic acid since this might trigger intolerances.

UCderm impression must not be mixed with other products.

The patient should avoid using any makeup for at least 12 hours after treatment as well as prolonged exposure to the sun, UV rays or temperatures below 0°C for at least 14 days after treatment. In addition, any sauna or hammam sessions should be avoided for up to two weeks after the injection.

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. Therefore, UCderm impression must never get into contact with these substances or with medical instruments which has been treated with these substances.

**Application:** Prior to treatment, inform the patient about the indication, contraindications, possible intolerances and potential side effects of UCderm impression.

This product is indicated for injection into the upper dermis by a legally qualified, medical practitioner. For successful treatment it is essential that the practitioner use the injection technique described in the enclosed injection instruction. He/she has to be knowledgeable about the anatomy and physiology of the areas to be treated. The treatment must be carried out under appropriate aseptic conditions. The area to be treated has to be disinfected thoroughly.

UCderm impression may be injected into healthy, non-inflamed skin areas only.

UCderm impression must not be injected into skin areas that contain residues of hyaluronic acid implants by other manufacturers. The last injection should be 3-6 months ago.

UCderm impression must not be injected into skin areas previously treated with an implant that is not based on hyaluronic acid.

Use the supplied, single-use sterile needle. For optimum use, assemble the needle as shown in Figures 1, 2 and 3 above.

- Firmly hold the glass cylinder of the syringe and the Luer-lock between the thumb and the index finger of one hand. Carefully remove the protective cap from the syringe with the other hand.
- Screw the supplied needle onto the syringe as far as it will go. Ensure not to twist the Luer lock.
- Keep holding the Luer-lock and carefully remove the protective cap of the needle.

UCderm impression should be injected slowly and precisely. Make sure to place the injection sites regularly and to treat the skin areas evenly. The amount to be injected depends on the correction to be achieved and the area to be treated. The syringe has an indexing scale.

Do not inject into blood vessels.

If the needle is blocked, do not increase the pressure on the plunger. Instead, stop the injection immediately and use a new needle. Never bend a straight needle for injection. Never try to straighten a bent needle. Discard the bent needle and replace it. If a needle appears dull, discard the blunt needle and use a new needle.

Gently massage the treated area according to the enclosed injection instruction after the injection to distribute the product uniformly.

The patient should avoid using any makeup for at least 12 hours after treatment as well as prolonged exposure to the sun, UV rays or temperatures below 0°C for at least 14 days after treatment. In addition, any sauna or hammam sessions should be avoided for up to two weeks after the injection.

Consider that the product contains substances which may produce a positive result in anti-doping tests.

**Side effects:** The patient must be informed that they are possible side effects, which may occur immediately or may be delayed. These include, but are not limited to:

- Slight bleeding during the treatment which spontaneously subsides as soon as the injection is completed.
- Inflammatory reactions usually associated with injections such as redness, oedema, erythema, swelling, pain or hyperthermia sometimes accompanied with itching, burning or paraesthesia. These reactions may last for 3-5 days.
- Delayed injection-related inflammation such as persistent erythema or acne papule formation.
- Bleeding, haematomas or telangiectasia in the treated area.

Swelling effect in other areas through displacement of adjacent structures or product migration.

– Immediate or delayed hypersensitivity or allergy (urticaria, in single cases anaphylaxis or Quincke’s edema) to one of the product’s components, especially to sodium hyaluronate and lidocaine hydrochloride.

- Infections, abscess, biofilm formation or recurrent herpetic lesions.
- Non-inflammatory and inflammatory nodule formation, induration or asymmetry.
- Skin discolouring at the injection site:
  - Due to superficial injection leading to bluish discoloration known as Tyndall Effect.
  - Due to ischemia leading to redness if vein is occluded or whiteness if artery is occluded.
  - Hyperpigmentation particularly in patients with Fitzpatrick skin type IV to VI.
- Non-inflammatory and inflammatory granuloma formation (foreign body reactions).
- Dermal fibrosis, ulcer, scarring or cellulitis.
- Lack of efficacy or reduced effect.

Cases of necrosis, muscle degeneration and osteitis after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

As a matter of general principle, injection of a medical device is associated with a risk of infection. Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe appropriate treatment.

**Warnings:**

- Check the integrity of the package and the expiry date prior to use. Do not use the product if the package has already been opened or is visibly damaged or expiration date has been exceeded. UCderm impression has a shelf life of 3 years after the date of manufacture.
- Do not re-sterilise the product.
- Do not re-use the product.
- Used and partially used syringes and injection needles should be handled and disposed in accordance with facility medical practices and applicable regulations. Do not remove the needle from the syringe or put the protective cap back on the needle.

**Storage:** Store UCderm impression between 2°C and 25°C and protect it from freezing and direct sunlight. Avoid mechanical shocks.

**Manufacturer:** UniCare GmbH Mittelstrasse 7 12529 Schoenefeld, Germany Website: www.unicare.berlin Phone:+49 (0) 30 63415601

**Item number (IFU):** UDBZ00116001 **Revision (IFU):** 3 (2019-07) **UDI-DI (UCderm impression):** 04260632410004 Please note: The summary of safety and clinical performance of UCderm impression is deposited at Eudamed and can be found via the UDI-DI.

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