



<div><div>U</div>derm lift</div> <div>1 x 1 ml</div>

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EN

Description: UCderm lift is a sterile, isotonic, pyrogen-free, colourless and transparent gel made of highly purified, crosslinked hyaluronic acid of non-animal origin in a physiological phosphat buffer. UCderm lift contains 0.3% lidocaine hydrochloride.

Content: Each package contains one single-use glass syringe, pre-filled with 1 ml UCderm lift and sterilised by moist heat. In addition, there is one single-use 27G sterile needle to be used only for injecting UCderm lift, an instruction leaflet, a patient card and two labels in order to ensure traceability in it.

Composition: UCderm lift contains:

- 26 mg/ml crosslinked Hyaluronic acid

- 3 mg/ml Lidocaine hydrochloride
- Phosphate buffer

The amount of gel in each syringe is 1 ml.

Indications/Intended use: UCderm lift is a sterile, injectable, absorbable implant for augmenting the volume of facial soft tissues. This implant is indicated to lift different types of deep facial wrinkles, to restore facial volume, to refine scars and to balance anatomical or injury-related asymmetries. Moreover, it is also intended for correction of facial fat loss caused by illness, medication or predisposition. The product is intended to be used for injection into the mid to deep dermis.

The lidocaine contained in the product is meant to reduce the patient’s local pain during and for a short time after the injection.

Contraindications: Do not apply UCderm lift in the following cases:

- In youths under 18 years of age.
- During pregnancy and breast-feeding.
- In case of incipient or existing influenza infections.
- In case of acute systemic infections.
- In case of skin inflammations (e.g. acne, herpes, psoriasis)
- In case of a known hypersensitivity to sodium hyaluronate or lidocaine hydrochloride or amide-based local anaesthetics.
- In case of known autoimmune diseases.
- In patients with blood clotting disorders or patients on anti-coagulation medication (e.g. aspirin, non-steroidal anti-inflammatory drugs).
- In case of acute articular rheumatism with cardiac localization.
- In patients with a tendency to hypertrophic scarring
- In case of untreated epilepsy.
- In patients who suffer from porphyria.
- During treatment with drugs tending to slow down or inhibit the liver metabolism (e.g. cimetidine, beta-blockers).
- In case of cardiac arrhythmia.

Precautions for use: UCderm lift must not be injected in the periorbital area (eyelids, eye shadows, crow’s feet) and glabellar region (forehead).

When treating facial areas with limited soft tissue support or soft tissue coverage, special care should be taken to avoid the formation of palpable lumps, discoloration or granulomas.

UCderm lift must not be injected into tendons, muscles, bones, joints or eyes.

UCderm lift is not intended for use in breast augmentation or breast reconstruction.

UCderm lift must not be injected into blood vessels. Introduction of product into the vasculature may lead to embolization or arterial ischemia or infarction. Therefore, the product should always be injected slowly and at the lowest possible pressure.

Do not overcorrect.

UCderm lift must not be used simultaneously with other aesthetic medicine techniques such as laser treatments, chemical peelings or dermabrasion. If inflamed skin areas occur after the treatment with these techniques, injection into those areas is not recommended. The treated skin area has to heal completely first.

UCderm lift must not be injected into skin areas which had been subjected to a treatment with permanent implants, implants of animal origin or implants containing other substances than hyaluronic acid since this might trigger intolerances.

UCderm lift must not be mixed with other products.

The patient should avoid using any makeup for at least 12 hours after treatment as well as prolonged exposure to the sun, UV rays or temperatures below 0°C for at least 14 days after treatment. In addition, any sauna or hammam sessions should be avoided for up to two weeks after the injection.

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride.

Therefore, UCderm lift must never get into contact with these substances or with medical instruments which has been treated with these substances.

Application: Prior to treatment, inform the patient about the indication, contraindications, possible intolerances and potential side effects of UCderm lift.

This product is indicated for injection into the mid to deep dermis by a legally qualified, medical practitioner. For successful treatment it is essential that the practitioner has been trained in injection technique for soft tissue augmentation. He/she has to be knowledgeable about the anatomy and physiology of the areas to be treated. The treatment must be carried out under appropriate aseptic conditions. The area to be treated has to be disinfected thoroughly.

UCderm lift may be injected into healthy, non-inflamed skin areas only.

UCderm lift must not be injected into skin areas that contain residues of hyaluronic acid implants by other manufacturers. The last injection should be 3-6 months ago.

UCderm lift must not be injected into skin areas previously treated with an implant that is not based on hyaluronic acid.

Use the supplied, single-use sterile needle. For optimum use, assemble the needle as shown in Figures 1, 2 and 3 above.

1. Firmly hold the glass cylinder of the syringe and the Luer-lock between the thumb and the index finger of one hand. Carefully remove the protective cap from the syringe with the other hand.

2. Screw the supplied needle onto the syringe as far as it will go. Ensure not to twist the Luer lock.

3. Keep holding the Luer-lock and carefully remove the protective cap of the needle.

UCderm lift should be injected slowly and evenly. The amount to be injected depends on the correction to be achieved and the area to be treated. Do not overcorrect. The syringe has an indexing scale.

If the needle is blocked, do not increase the pressure on the plunger. Instead, stop the injection immediately and use a new needle.

Never bend a straight needle for injection. Never try to straighten a bent needle. Discard the bent needle and replace it.

If a needle appears dull, discard the blunt needle and use a new needle.

Gently massage the treated area after the injection to distribute the product uniformly.

The patient should avoid using any makeup for at least 12 hours after treatment as well as prolonged exposure to the sun, UV rays or temperatures below 0°C for at least 14 days after treatment. In addition, any sauna or hammam sessions should be avoided for up to two weeks after the injection.

Consider that the product contains substances which may produce a positive result in anti-doping tests.

Side effects: The patient must be informed that they are possible side effects, which may occur immediately or may be delayed. These include, but are not limited to:

- Slight bleeding during the treatment which spontaneously subsides as soon as the injection is completed.
- Inflammatory reactions usually associated with injections such as redness, oedema, erythema, swelling, pain or hyperthermia sometimes accompanied with itching, burning or paraesthesia. These reactions may last for 3-5 days.
- Delayed injection-related inflammation such as persistent erythema or acne papule formation.
- Bleeding, haematomas or telangiectasia in the treated area.
- Swelling effect in other areas through displacement of adjacent structures or product migration.

- Immediate or delayed hypersensitivity or allergy (urticaria, in single cases anaphylaxis or Quincke’s edema) to one of the product’s components, especially to sodium hyaluronate and lidocaine hydrochloride.
- Infections, abscess, biofilm formation or recurrent herpetic lesions.
- Non-inflammatory and inflammatory nodule formation, induration or asymmetry.
- Skin discolouring at the injection site:
 - Due to superficial injection leading to bluish discoloration known as Tyndall Effect.
 - Due to ischemia leading to redness if vein is occluded or whiteness if artery is occluded.
- Hyperpigmentation particularly in patients with Fitzpatrick skin type IV to VI.

- Non-inflammatory and inflammatory granuloma formation (foreign body reactions).
- Dermal fibrosis, ulcer, scarring or cellulitis.
- Lack of efficacy or reduced effect.
- Cases of necrosis, muscle degeneration and osteitis after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

As a matter of general principle, injection of a medical device is associated with a risk of infection.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe appropriate treatment.

Warnings:

- Check the integrity of the package and the expiry date prior to use. Do not use the product if the package has already been opened or is visibly damaged or expiration date has been exceeded. UCderm lift has a shelf life of 3 years after the date of manufacture.
- Do not re-sterilise the product.
- Do not re-use the product.
- Used and partially used syringes and injection needles should be handled and disposed in accordance with facility medical practices and applicable regulations. Do not remove the needle from the syringe or put the protective cap back on the needle.

Storage: Store UCderm lift between 2°C and 25°C and protect it from freezing and direct sunlight. Avoid mechanical shocks.

Manufacturer:
UniCare GmbH
Mittelstrasse 7
12529 Schoenefeld, Germany
Website: www.unicare.berlin
Phone: +49 (0) 30 63415601

Item number (IFU): UDBZ00116004
Revision (IFU): 3 (2019-07)

UDI-DI (UCderm lift): 04260632410035

Please note: The summary of safety and clinical performance of UCderm lift is deposited at Eudamed and can be found via the UDI-DI.

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